

TECHNICAL DATA SHEET

Distributed by:

AMA Medical Products

amamedicalproducts.com.au

sales@amamedicalproducts.com.au 1800 626 292


Description

7000 Impervious Film Gown, Thumb hooks, Open back, Blue, Regular

7001 Impervious Film Gown, Thumb hooks, Open back, Blue, X-Large

7002 Impervious Film Gown, Thumb hooks, Open back, Yellow, Regular

7003 Impervious Film Gown, Thumb hooks, Open back, Yellow, X-Large



Dimensions

	Regular	X-Large
Overall Length	109 cm	124 cm
Body Width	93 cm	111 cm
Sleeves End-To-End	191 cm	211 cm
Neck Opening	24 cm	27 cm

Properties

Resistance of Materials to Penetration By Synthetic Blood – ASTM1670	Time to Penetration >60 mins	Pass
Resistance of Materials to Penetration By Blood-Borne Pathogens – ASTM1671	Time to Penetration >60 mins	Pass

Main materials

- Impervious polyethylene film.
- Does not contain Natural Rubber Latex. Does not contain DEHP.

Packaging

Shipping case of 75 units made up of:

- 15 units per dispenser box
- 5 boxes per case

Sterilization

These products are non-sterile.

Manufacturing

Made in Thailand.

The manufacturing site has ISO 13485 certified quality system.

Storage

Store in a dry and cool place, away from intense sources of heat. Keep the gowns as much as practicably possible in their dispenser box. Keep dispenser boxes as much as practicably possible in their shipper box.

Shelf life

Non-specified.

Indication

To minimise the direct transmission of infective agents between the healthcare staff and the patient. Single use product.

Potential for soaking and HEAVY FLUID CONTACT in front.

Contraindication

None in particular.