


Endotracheal Tube Cuffed

PRODUCT DATA SHEET



DESCRIPTION

The m|devices Endotracheal Tube - Cuffed are used to assist patients to perform respiration while connected to the respiratory system during a surgical procedure or who lose active breath ability. These ET tubes are available with an inflatable balloon (or cuff) to secure the tube in place.

FEATURES AND BENEFITS

- Transparent, soft and smooth
- Murphy eye design provides an alternate gas passage
- Radio opaque line for X-Ray visualisation
- High volume, positive cuff pressure
- One way valve for inflation
- Kink resistant
- Large size range
- Bold gradient markings

PRODUCT DATA

- Pilot Colour: Blue
- Latex: No
- DEHP: No
- Sterile: EO
- Materials: PVC
- Shelf life: 3 Years
- ARTG Number: 210994

RE-ORDER	DESCRIPTION	TUBE LENGTH	CARTON/QTY
AN021000	Endotracheal Tube - Cuffed, 3.00mm	190mm	10
AN021001	Endotracheal Tube - Cuffed, 3.5mm	210mm	10
AN021002	Endotracheal Tube - Cuffed, 4.00mm	230mm	10
AN021003	Endotracheal Tube - Cuffed, 4.5mm	240mm	10
AN021004	Endotracheal Tube - Cuffed, 5.0mm	260mm	10
AN021005	Endotracheal Tube - Cuffed, 5.5mm	270mm	10
AN021006	Endotracheal Tube - Cuffed, 6.0mm	290mm	10
AN021007	Endotracheal Tube - Cuffed, 6.5mm	310mm	10
AN021008	Endotracheal Tube - Cuffed, 7.0mm	320mm	10
AN021009	Endotracheal Tube - Cuffed, 7.5mm	330mm	10
AN021010	Endotracheal Tube - Cuffed, 8.0mm	340mm	10
AN021011	Endotracheal Tube - Cuffed, 8.5mm	350mm	10
AN021012	Endotracheal Tube - Cuffed, 9.0mm	350mm	10
AN021013	Endotracheal Tube - Cuffed, 9.5mm	350mm	10
AN021014	Endotracheal Tube - Cuffed, 10.0mm	350mm	10