



AQ Medical 280 Vaccine Refrigerator

The AQ Medical 280 vaccine refrigerator is a robust and cost effective unit, purpose-built for vaccine storage following the Australian Government Department of Health's 'Strive for 5' vaccine storage guidelines.

It was designed with low noise output and consistent temperature readings in mind, and maintains a stable temperature between the recommended +2°C to +8°C. The unit holds safety features such as visual and audible alarms in the event that temperatures fluctuate outside the set range.

All AQ Medical refrigerators are accredited by the Quality Care Pharmacy Program (QCPP), a quality assurance program that was developed by the Pharmacy Guild of Australia in consultation with the Pharmaceutical Society of Australia.

AQ Medical units are also backed by the extensive Avem Quirks service network. With qualified technicians and warehousing facilities across Australia, we can deliver responsive service and exceptional support on a national level.

The AQ Medical 280 vaccine refrigerator is registered with GEMS and has a 8 Star Rating, one of the most energy efficient vaccine refrigerators on the Australian Market. Registration Number: ARC1538



Distributed by:
AMA Medical Products
amamedicalproducts.com.au
sales@amamedicalproducts.com.au
1800 626 292



AQ Medical 280 Vaccine Refrigerator

Key Features

- QCPP cold chain accredited
- Digital temperature control system preconfigured to hold temperature from +2°C to +8° C in line with national 'Strive for 5' storage guidelines
- Audible and visual auto resetting alarms
- Triple glazed glass
- Adjustable shelves
- Lockable door
- Low energy consumption
- Low noise output
- Dedicated technical service line

Optional Extras

- Uninterruptible Power Supply unit with back up battery
- AQ Box cloud temperature monitoring system
- Comprehensive programmed maintenance

Product Specifications

Dimensions	Height 1575 mm
	Width 584 mm
	Depth 584 mm
Gross Weight	79 kg
Gross Capacity	280 L
Temperature Range	+2°C to +8°C
Number of Shelves	6
Power Consumption/24 hrs.	2.6 kWh
Refrigerant	R600a
Warranty (parts and labour)	24 months



Distributed by:
AMA Medical Products
amamedicalproducts.com.au
sales@amamedicalproducts.com.au
1800 626 292

