



Smart way to a healthy life



connects with
HealthVault
Microsoft



Continua™
CERTIFIED

CE 0123

A&D
A&D Medical

Technology You Can Trust

<http://www.aandd.jp>
<http://www.wellnessconnected.jp/english/>

Digital Blood Pressure Monitor UA-651BLE



Precision Health Scale UC-352BLE



- Clinically validated for accuracy (ESH2010)
- Bluetooth Ver.4.0LE
- Bluetooth communication (Windows, iOS, Android)
- Internal memory of 30 data records
- Irregular HeartBeat (IHB) indicator - LCD only
- Extra large 3-line display for easier reading of results
- One-touch measurement
- Latex-free and metal-free, slimFit Cuff
- Compact design for portability

- Bluetooth Ver.4.0LE
- Bluetooth communication (Windows, iOS, Android)
- Durable glass plate scale
- kg/lb selection
- Internal memory of 90 data records
- 0 to 100kg : 100g resolution / 100 to 200kg : 200g resolution
- 0 to 220.0lb : 0.2lb resolution / 220.0 to 450.0lb : 0.5lb resolution
- Simply step on start weighing
- Beeps on completion of weighing
- Low profile stylish design

Specifications

UA-651BLE

Measurement method	Oscillometric measurement
Display type	Digital LCD, Pressure/Pulse displayed simultaneously
Measurement range	Pressure: 0 - 299 mmHg Systolic pressure: 60 - 279 mmHg Diastolic pressure: 40 - 200 mmHg Pulse: 40 - 180 beats / minute
Measurement accuracy	Pressure: ± 3 mmHg Pulse: $\pm 5\%$
Power supply	4 x 1.5V batteries (LR6, R6P or AA) or AC adapter (TB-233) (Not included)
Number of measurements	Approx. 700 times LR6 (alkaline batteries) Approx. 200 times R6P (manganese batteries) With pressure value of 180 mmHg at room temperature of 23 °C
Memory storage	30
Cuff	SlimFit Cuff with plastic D-ring / Arm circumference : 22-32cm
Operating conditions	+10°C to +40°C / 15%RH to 85%RH/800 hPa to 1060 hPa
Transport/Storage conditions	-20°C to +60°C / 10%RH to 95%RH
Wireless communication	Bluetooth Ver.4.0LE BLP
Communication specifications	*Continua Design Guideline 2013 Version 4.1 (ITU-T H.810 (12/2013) Part 0, Part 1)
Dimensions	Approx. 96 [W] x 68 [H] x 130 [D] mm
Weight	Approx. 250 g, excluding batteries
Options	SlimFit Large Cuff (31-45cm) : CUF-F-LA Wide Range Cuff (22-42cm) : CUF-I AC Adapter (C/BF plug) : TB-233C/BF

UC-352BLE

Maximum capacity	200 kg / 450 lb
Display type	Digital LCD, Weight displayed after step off
Minimum display	0.1 kg (0 to 100 kg) / 0.2 kg (100 to 200 kg) 0.2 lb (0 to 220 lb) / 0.5 lb (220 to 450 lb)
Setting items	Units: kg or lb
Power source	4 x 1.5V alkaline batteries (LR03 or AAA)
Battery life	Approximately 1 year (two measurements per day when alkaline batteries are used)
Memory storage	90
Operating conditions	+10°C to +40°C / 30%RH to 85%RH
Transport/Storage conditions	-10°C to +60°C / 30%RH to 85%RH
Wireless communication	Bluetooth Ver.4.0LE WSP
Communication specifications	*Continua Design Guideline 2014 Version 5.0
Dimensions	Approx. 300 [W] x 300[D] x 31[H] mm
Weight	Approx. 2.1 kg (4.7lb) including batteries

*Continua certification is instituted and promoted by PCHA to certify compliance to Continua Design Guidelines. The guidelines were developed to realize an end-to-end interconnected service for healthcare products.

The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by A&D is under license. Other trademarks and trade names are those of their respective owners. WellnessConnected™ is trademark of A&D Company, Limited or A&D Medical. Apple name and logo are registered trademarks of Apple, Inc. Google, Android, Google Logo and Google Play are registered trademarks of Google, Inc. Facebook name and logo are trademarks of Facebook Inc. Twitter name and logo are trademarks of Twitter Inc. Microsoft HealthVault icon is a registered logo of Microsoft Corporation in the United States and/or other countries.


A&D Medical is preparing a range of healthcare devices compatible with A&D Connect app.




A&D Connect

Easy-to-use app, a solution for better accuracy


Monitor your blood pressure and track your weight and BMI trends – all with one app from a trusted world-wide leader in connected health. A&D Connect guides you through the easy device set-up process and provides at-a-glance results, goals, and trends for you and sharing data with your family and doctors.




Dashboard



List



Graph



Sharing



Available on iOS 6.0 and above
iPhone 4S or newer, iPad Air, iPad mini



Requires Android OS version 4.4 and above
For compatible smartphones, please check WellnessConnected™ website



A&D Company, Limited
3-23-14 Higashi-Ikebukuro, Toshima-ku, Tokyo 170-0013 Japan
Telephone: [81](3)5391-6132 Fax: [81](3)5391-6148
1-243, Asahi, Kitamoto-shi, Saitama, 364-8585 Japan
Telephone: [81](0485)93-1111 Fax: [81](0485)93-1119
<http://www.aandj.jp>

A&D INSTRUMENTS LIMITED
Unit 24/26 Blacklands Way Abingdon Business Park,
Abingdon, Oxon OX14 1DY UNITED KINGDOM
Telephone: [44](1235) 550420 Fax: [44](1235) 550485
(Authorized Representative Established in the European Community)

A&D INSTRUMENTS INDIA PRIVATE LIMITED
509 Udyog Vihar Phase V Gurgaon-122 016, Haryana, India
Telephone : [91](124)471-5555 Fax: [91](124)471-5599