

A&D Medical ... *Clearly a Better Value*

PRECISION HEALTH SCALE UC-355PBT-Ci



The UC-355PBT-Ci Precision Health Scale has a 550lb capacity, supporting a large range of users. This scale can display your weight in pounds or kilograms, allowing you to select whichever you are comfortable with. Additionally, the scale has a wider platform to support those with a wide base. It is also equipped with audible buzzers and visual indicators to let you know when to get on and off. The profile is less than 2" tall even when the carpet tray (included) is utilized. You don't have to push any buttons to turn on this scale, with A&D's step-on technology. Just step on and let the scale guide you with on-screen prompts.

The PBT-Ci family of products boasts updated Bluetooth technology (v2.1) which allows for improved connectivity through secure simple pairing (SSP). This device can operate in various profiles (e.g. Serial Port Profile (SPP), Health Device Profile (HDP) and gives you the choice between master or slave mode operation. System integrators are even provided enhanced data with measurement status (e.g. low battery, etc).

For questions about A&D products, call us toll free at 1-888-726-9966 and a specially trained representative will assist you, or visit us at www.andmedical.com.



- Class 1 Bluetooth
- High Capacity - 550 lb (250kg)
- Low Profile - 1.6"
- Step-on technology
- Large Display - 35mm
- Includes Scale tray
- 200 measurement memory

PRECISION HEALTH SCALE UC-355PBT-Ci

PART NUMBER	UPC NUMBER	DESCRIPTION
UC-355PBT-Ci	5-060006-580932	550lb Precision Health Scale



Technical Specifications

Measuring method	Electric resistance
Digital Display	Liquid crystal display, 35 mm
Maximum capacity	550 lbs / 250 kg
Resolution	0.1 kg between 0 - 100 kg / 0.2 lb between 0 - 220 lb 0.2 kg between 100 - 250 kg / 0.5 lb between 220 - 550 lb
Power Source	Four type AA batteries (1.5 volt)
Dimensions	14.2" (d) x 16.4" (w) x 1.5" (h) / 360 mm (d) x 415 mm (w) x 39 mm (h)
Weight	Approximately 10.0 lbs (4.6 kg), including the batteries
Accessories	Instruction manual, batteries (AA x 4), scale tray
Battery Life	Approximately 1 year with two measurements per day
Setting Items (Units)	kg or lb
Communication	Bluetooth Ver.2.1 class 1 SSP HDP compatible

Packaging and Shipping Specifications

Selling Unit (Box)	Dimensions (L x W x H)	Weight	Individual Item Cube (cubic feet)
UC-355PBT-Ci	22 x 4 x 20.5 in	12.6 lb	1.04
Unit Case (3 units)	Dimensions (W x D x H)	Weight	Case Cube (cubic feet)
UC-355PBT-Ci	23 x 20 x 10 in	42 lb	2.67

Ask your A&D Medical Sales Representative about other A&D Medical products.

- Blood Pressure Monitors
- Digital Thermometers
- Personal Scales
- Activity Monitors



1756 Automation Parkway
San Jose, CA 95131
www.andmedical.com

© 2013 A&D Medical. All rights reserved.

For questions about A&D products, call us toll free at 1-888-726-9966 and a specially trained representative will assist you, or visit us at www.andmedical.com.