



DIAGNOSTIC POCKETSCOPE 2.5v AA SET

STOCK CODE: WEL92820

❖ PocketScope Ophthalmoscope

- Lightweight, compact diagnostic instruments deliver high performance in a slightly smaller package
- Halogen light for true tissue colour and consistent, long-lasting illumination
- 12 aperture/filter combinations for greater versatility: micro, small, and large spot sizes, cobalt blue filter for corneal exams, fixation target, and slit aperture, plus re-free filter & unfiltered halogen
- 8 focusing lenses for better resolution
- PocketScope Otoscope with built-in throat Illuminator
- Small, convenient, and easy to use with bright, white halogen light and a built-in throat illuminator
- Black hard case
- Sample ear spec
- Powered by AA batteries
- Metal handle casing

PocketScope Ophthalmoscope Accessories

- 2.5 v Halogen Lamp Stock Code: WEL03900
- Black Hard Case Stock Code: WEL92820
- Black Soft Case Stock Code: WEL92821

PocketScope Otoscope / Throat Illuminator Accessories

- 2.5 V Halogen Lamp Stock Code: WEL03400
- Insufflation Bulb with Tip for Diagnostic Otoscope Stock Code: WEL21504