



Welch Allyn  
CP 150™ Electrocardiograph  
with optional Spirometry

Simple. Fast. Connected.



# Keep your practice running with the Welch Allyn CP 150 Resting ECG

Simple, streamlined ECG testing with a range of advanced features to help improve your workflow.

The CP 150 combines the best traditional-box ECG features like full-size printouts and clinical decision support with MEANS interpretation with a simple touch-screen interface, fast one-button operation, and flexible connectivity options to help you keep your staff productive and your tests accurate. In fact, 81% of clinicians prefer the CP 150 over the competition.\*

\*Source: Walker Information, Inc.




## Simple



The 7-inch, touch-screen display with on-screen prompts guides staff through a test, to help enhance accuracy, minimize training and save time.



Avoid repeat ECGs with the lead quality display that is designed to help technicians quickly and easily identify lead connection problems before completing a test.



of CP 150 users say the touchscreen display and simple user interface has helped them improve their practice efficiency.


Source: Survey by Walker Information, Inc. of clinicians who have purchased the CP 150.

## Fast




### Instant-On Feature

Powers the device up quickly, and eliminates the typical startup delay to take an ECG.



### 12-Lead ECG Preview

Large display provides three, six, or 12-lead ECG preview options.




### ECG Print Preview

Save time and reduce paper costs by verifying the ECG test quality prior to printout.



## Spirometry Option

The CP 150 lets you perform two tests with one device, to help make detection at the point of care easier and more efficient.

- > Perform FVC and pre- and post-BD testing with automatic interpretation and reversibility analysis
- > Customize report formats and screen displays to fit your needs and workflow
- > Performs spirometry tests typically required for health and safety purposes including ATS and NIOSH
- > Automatic lung age calculation to help with smoking cessation programs
- > Integrated ambient sensors that eliminate the need for the user to find and enter values for temperature, pressure and humidity



The CP 150 spirometry option features four incentive screens to help educate and coach your patients:



# Connectivity Options that mean you'll never have to scan and shred an ECG test again.



The CP 150 features a range of connectivity options to match your workflow



## Simple, Streamlined EMR Integration

- Save ECG tests in PDF format
- Export tests to USB and drag-and-drop file into folder in EMR
- Use the integrated ethernet connection and export PDF directly to a folder on your network



## Bring your cardiopulmonary data together with the Welch Allyn CardioPerfect® Workstation

- Store, access, review, compare and edit diagnostic test data from any computer on your network, in any exam room, or outside the office
- Ability to interface with more than 40 Electronic Health Record systems

## Take your CP 150 on the go

- Attractive, ergonomic mobile cart features optional hospital grade wheels and patient cable arm and shelf
- Long-lasting lithium ion battery gives you more on-time, letting you transport your device directly to the patient



## Ordering Information

### Device Configurations:

CP150A-1ENB	CP 150 Interpretive Electrocardiograph
CP150-1ENB	CP 150 Non-Interpretive Electrograph
CP150AS-1ENB	CP 150 Interpretive Electrocardiograph with Spirometry
CP150S-1ENB	CP 150 Non-Interpretive Electrocardiograph with Spirometry

### Accessories:

715006	ECG Multifunction Electrode Adaptor
45008-0000	Resting Tab Electrodes (box of 1000)
105353	CP 150 Recording Paper (case, 1000 sheets)
719653	10-lead Patient Cable, AHA, Banana, CP 150
105341	CP 150 Office Cart
105342	CP 150 Hospital Cart
105343	CP 150 Cable Arm & Shelf Cart

### Spirometry Accessories:

105660	Spirometry Upgrade for CP 150 (works with any standard CP150A-1ENB or CP150-1ENB configuration)
720705	Disposable Flow Transducers, 25 pk, CP 150
720706	Disposable Flow Transducers, 100 pk, CP 150 / Spiro

### Partners in Care<sup>SM</sup> Service Options:

S1-CP150	CP 150 Comprehensive Partner Program
S-1-CP150-2	CP 150 Comprehensive Partner Program, 2 YR
S1-CP150S	CP 150 Spiro Comprehensive Partner Program
S1-CP150S-2	CP 150 Spiro Comprehensive Partner Program, 2 YR

### Warranty Information:

CP 150 Electrocardiograph	—3 years
Patient Cable/Battery	—90 days
Spirometer	—1 year

Contact your Welch Allyn representative or visit [www.welchallyn.com/cp150](http://www.welchallyn.com/cp150) today for more information.

Welch Allyn Corporate Headquarters  
4341 State Street Road, P.O. Box 220  
Skaneateles Falls, NY 13153-0220 USA  
(p) 800.535.6663 (f) 315.685.3361

[WWW.WELCHALLYN.COM](http://WWW.WELCHALLYN.COM)

**WelchAllyn**<sup>®</sup>