

Introducing the **NEW** Welch Allyn SureColor™ LED Lamps



The new Welch Allyn SureColor™ LED lamp provides a greater than **80% cost saving** over the life of the lamp versus our equivalent halogen lamp, plus:

3X on-time for rechargeable battery handles

Lower heat build up for increased clinician and patient comfort

Consistent color output through dimming of light rheostat

Patented circuitry and color saturation for specific clinical application

7-year warranty

Average 30X longer light hours for approximately 3X the price

WelchAllyn®

Advancing Frontline Care™

Welch Allyn SureColor™ LED lamps offer:

SureColor™ LED Technology enabling reduction of light intensity to alleviate glare and enhance patient comfort while preserving the color and clarity of the view.

High Energy Efficiency resulting in 3X the examination time with each charge cycle of a Welch Allyn rechargeable handle.

A Simple Solution presented with four lamp upgrade kits which are compatible with our most popular physical assessment instruments and can be installed in seconds with a simple exchange of the halogen lamp.

LED Lamp Part#		Device Part#	Description
Traditional Diagnostic Oscopes			
03100-LED Halogen: 03100-U 		20000	3.5V Oscope/Throat Illuminator
		20200	3.5V Pneumatic Oscope
		21700	3.5V Operating Oscope
		25020	3.5V Diagnostic Oscope (with Specula)
MacroView™ Diagnostic Oscopes			
06500-LED Halogen: 06500-U 		23810	3.5V MacroView™ Oscope
		23820	3.5V MacroView™ Oscope/Throat Illuminator
		23920	Digital MacroView Oscope
Traditional Ophthalmoscopes			
04900-LED Halogen: 04900-U 		11720	3.5V Coaxial Ophthalmoscope
		11730	AutoStep® Coaxial Ophthalmoscope
		11735	3.5V Prestige™ Coaxial-Plus Ophthalmoscope
PanOptic™ Ophthalmoscopes			
03800-LED Halogen: 03800-U 		11820	PanOptic™ Ophthalmoscope/Blue Filter/Lens

*Prices exclude 10% GST. LED Lamp Upgrade Kit 3.5V Diag [10/box] also available.

For more information contact Welch Allyn your Welch Allyn Representative on 1800 650 083
www.welchallyn.com/ledlamps