

LGex 3410

MediLine



Total gross capacity	310 litre
Ext. dimensions in mm (w/d/h)	600/615/1840
Int. dimensions in mm (w/d/h)	420/400/1587
Energy consumption in 24 h	1,309 kWh
Climate rating ¹	SN-ST
Cooling system	static
Defrosting	manual
Temperature range	-9°C to -30°C
Gastronorm	-
Housing material/colour	steel/white
Door/lid material	Steel
Liner Material	moulded polystyrol white
Type of control	electronic controller
Temperature display	external digital
Interior light	-
Fault: warning signal	visual and audible
Interface	RS 485
volt-free alarm contact	yes
Shelves	8
Adjustable Shelving	-
Number of drawers	8
Number of baskets	-
Maximum number of baskets	8
Shelf material	evaporator plates
Shelf loading	24 kg
Height of adjustable legs	-
Castors/type of castors	-
Door handle	Ergonomic door handle
Foot pedal opener	-
Self-closing door	yes
Insulation	70 - 70 mm
Defrost water drain	-
Refrigerant	R 600a
Noise level	45 dB(A)
Gross/net weight	92/87 kg
EAN-Nr.	9005382171059
Connection Rating	1.5 A / 280 W

¹ Climate Ratings: **SN** +10°C to +32°C **N** +16°C to +32°C **ST** +16°C to +38°C **T** +16°C to +43°C