

Welch Allyn Green Series™ 300 General Exam Light



Touchless On/Off



Easy Power Connector



Easy Assembly



1 LED

Whiter, Brighter, and Greener

The GS 300 General Exam Light features:

- A single LED that provides whiter, brighter light than halogen
- Less energy consumption for reduced cost of ownership
- No bulb replacement for cost savings
- A perfect combination of intense light output and cool operation, all delivered in a compact design
- An affordable price and rugged performance
- Minimal heat dissipation
- Touchless on/off
- 5-year limited warranty

WelchAllyn®
Advancing Frontline Care™

Green Series Medical Exam Lights



Welch Allyn Green Series™ 300 General Exam Light

ORDERING INFORMATION

- 44400 GS 300 General Exam Light with Mobile Stand
- 44410 GS 300 General Exam Light with Table/Wall Mount
- 44215 12" Wall Mount Extension
- 48960 Heavy Duty/Tall Mobile Stand

SPECIFICATIONS


# of LEDs	1 White LED
Color Temperature	5,500° K
Light Intensity @ Typical Working Distance	20,000 lux (1,858 footcandles)
LED Life	50,000 hours
Working Distance	12 in (30.5 cm) to 20 in (50 cm)
Spot Size at Working Distance	5 in (13 cm) to 6 in (15 cm)
Electrical	Input: 100-240 V~ 50 Hz-60 Hz, 0.35 A max
Power Cord	#18 AWG, 3-wire grounded cord set, 8' (2.5 m) long
Physical Dimensions	Weight: Total Table/Wall Mount: 5.65 lbs (2.56 kg) Total Mobile Stand: 12.75 lbs (5.78 kg) Control Box: Height: 10.83 in (27.5 cm) Width: 3.94 in (10.0 cm) Depth: 3.15 in (8.0 cm) Luminaire: Height: 7.48 in (19 cm) Width: 3.54 in (9 cm) Depth: 3.15 in (8 cm) Mobile Stand: Height: 24.02 in (61 cm) Diameter: 20.87 in (53 cm) Total Wall Mount: 33.27 in (84.5 cm) Total Mobile Stand: 58.66 in (149.0 cm) Table/Wall and Mobile Stand to Luminaire Center: 31.1 in (79 cm)
Mounting Options	Table/Wall Mount Mobile Stand
Control Box	Touchless On/Off Control
Accessory Options	NA
Warranty	5-year limited



IDEAL FOR USE IN THE FOLLOWING ENVIRONMENTS:

- PRIMARY CARE FACILITIES
- EMERGENCY ROOMS
- OUTPATIENT FACILITIES
- DERMATOLOGY
- URGENT CARE CENTERS

MOUNTING OPTIONS AVAILABLE:



To place an order, contact our Customer Service Department at 1.800.535.6663



WelchAllyn
Advancing Frontline Care™

