




# Dewars

Distributed by:  
 AMA Medical Products  
 amamedicalproducts.com.au  
 sales@amamedicalproducts.com.au  
 1800 626 292

## Liquid Nitrogen Container - Dewar



- More than 5 year's life
- Made of high-strength aluminium alloy
- Adopting advanced powder coating technology
- Each Dewar unit includes a leather (PU) Bag & Ladle (Pic. C)
- Optional extra's include: Polar gloves (Pic. A), Extractor (Pic. B)

Model		Unit	10L	20L	30L	
Capacity		L	10	20	30	
Weight Empty		kg	6.0	9.5	12.9	
Caliber		mm	50	50	50	
Outer Diameter		mm	300	394	462	
Height		mm	552	672	706	
Static evaporation time rate		L/d	0.12	0.12	0.12	
Static Holding Time		days	86	168	254	
Canister	Outer Diameter	mm	38	38	38	
	Height	mm	120	120/276	120/276	
	Canisters No. (pc)	pc	3-6	3-6	3-6	
	Canister No. (Standard) (pc)	pc	1	1	1	
	Straws holding capacity (single 120mm)	0.5ml	pc	792	792	792
		0.25ml	pc	1788	1788	1788
	Straws holding capacity (double 276mm)	0.5ml	pc	-	1284	1284
		0.25ml	pc	-	2832	2832

Distributed by:  
**AMA Medical Products**  
 amamedicalproducts.com.au  
 sales@amamedicalproducts.com.au  
 1800 626 292