

Wheelchair - Standard

Product Description

- Frame: Single cross brace, powder coating
- Front wheel: 8"/6", rubber tires
- Rear wheel: 24" spoke wheel, pneumatic tires
- Armrest: Flip back with desk length
- Footrest: Swing away footrest
- Front fork: Steel
- Footplate: Plastic



Size	Overall Length (cm)	Overall Width (cm)	Overall Height (cm)	Overall Width Folded	Seat Width (cm)	Seat Depth (cm)	Back Height (cm)	Floor to seat Height	Weight Capacity (kg)	Front Caster Spec.	Rear Wheel Spec.	G/N Weight (kg)	Carton Size (cm)
18'	103	63	91	26.5	45	40	44	50	110	8'	24'	22.5/20	82x27x43

Wheelchair - Patient mover

Product Description

- Frame: Single cross brace, powder coating
- Front wheel: 8", PVC tires
- Rear wheel: 12" PU tires
- Armrest: Flip back with desk length
- Footrest: Swing away footrest
- Front fork: Steel
- Footplate: Plastic



Size	Overall Length (cm)	Overall Width (cm)	Overall Height (cm)	Overall Width Folded	Seat Width (cm)	Seat Depth (cm)	Back Height (cm)	Floor to seat Height	Weight Capacity (kg)	Front Caster Spec.	Rear Wheel Spec.	G/N Weight (kg)	Carton Size (cm)
18'	97	63.5	92	27	46	41	43	49	110	8'	12'	16.4/14.4	72x25x93

Wheelchair - Bariatric

Product Description

- Frame: Single cross brace, powder coating
- Front wheel: 8"x2", PU tires
- Rear wheel: 24" spoke wheel, PU tires
- Armrest: Flip back with desk length
- Footrest: Swing away footrest
- Seat and backrest: Nylon with sponge
- Front fork: Steel
- Footplate: Plastic



Size	Overall Length (cm)	Overall Width (cm)	Overall Height (cm)	Overall Width Folded	Seat Width (cm)	Seat Depth (cm)	Back Height (cm)	Floor to seat Height	Weight Capacity (kg)	Front Caster Spec.	Rear Wheel Spec.	G/N Weight (kg)	Carton Size (cm)
22'	103	73	91	26.5	55	40	44	50	160	8"	24"	23/20.5	82x27x93