




## Athena Surgical Face Mask



The PrimeOn Athena Surgical Face Mask provides protection at a premium level for surgical tasks. Our anti-fog foam technology prevents condensation from escaping out of the top of the mask, stopping eyewear from fogging, and ensuring maximum visibility at all times. The PrimeOn Athena surgical face mask is resistant against fluid splashes and conforms to the Australian Standard. The tie design ensures comfort for different head shapes and sizes.



PrimeOn Athena with Ties	
TGA Registration	ARTG 290961
Application	Ties
Number of Fabric Layers	4
Anti-fog Foam	Yes
Latex Free	Yes
Colour	White
Fluid Resistance (mmHg)	160
Bacterial Filtration Efficiency	>99%
Differential Pressure (Delta Δ P), mm H <sub>2</sub> O/cm <sup>2</sup>	<5.0

**REORDER CODE**  
 FPM33UN

**LEVEL BARRIER CLASSIFICATION**  
 Level 3

**RECOMMENDED USE**  
 Where the risk of exposure to blood or bodily fluid splash are high

**PACKAGING**  
 50 masks per box  
 6 boxes per carton

**STANDARDS**  
 AS4381-2015, ASTM F2100

**MANUFACTURING ACCREDITATIONS**  
 ISO 9001, ISO 13485

A brand by



**Mun Australia**  
 Toll free: 1800 456 837  
 munglobal.com.au



Proud Supporter of  
 Birthing Kit Foundation (Australia)



Protection Always On