

# LKv 3910

Laboratory refrigerator with  
Comfort controller  
**MediLine**



Gross/net capacity	361/333 litre
Ext. dimensions in mm (w/d/h)	601/618/1840
Int. dimensions in mm (w/d/h)	440/435/1635
Energy consumption in 365 days	309 kWh
Ambient temperature range	+10 °C to +35 °C
Refrigerant	R 600a
Connection Rating	1.5 A
Frequency / voltage	50 Hz / 220-240 V~
Cooling system	forced-air
Defrosting	automatic
Temperature range	+3 °C to +16 °C
Housing material/colour	steel/white
Door material	Steel
Liner Material	moulded polystyrol white
Type of Control	electronic controller
Temperature display	external digital
Alarm in case of malfunction	optisch und akustisch
Interface	RS 485
volt-free alarm contact	yes
Adjustable Shelving	9
Shelf loading	45 kg
Lock	fitted
Self-closing door	yes
Door hinges	right, reversible
Insulation	40 mm
Gross/net weight	70/65 kg
EAN-Nr.	9005382237915

