

LCv 4010

MediLine



Gross/net capacity	
Fridge compartment	254/240 litre
Freezer compartment	107/105 litre
Ext. dimensions in mm (w/d/h)	
601/618/2003	
Int. dimensions in mm (w/d/h)	
Fridge compartment	440/441/1105
Freezer compartment	433/433/597
Energy consumption in 365 days	657 kWh
Ambient temperature range	+10°C to +35°C ¹
Refrigerant	R 600a
Connection Rating	1.5 A
Voltage	220-240 V~
Cooling system	forced-air
Defrosting	automatic
Temperature range	
Fridge compartment	+3°C to +16°C
Freezer compartment	-9°C to -30°C
Housing material/colour	steel/white
Door material	Steel
Liner Material	moulded polystyrol white
Type of Control	electronic controller
Temperature display	external digital/external digital
Fault: warning signal	visual and audible
Interface	RS 485
volt-free alarm contact	yes
Adjustable Shelving	4
Drawers	3
Shelves, freezer compartment	3
Compartment height in mm	2 × 156, 1 × 165
Lock	fitted
Self-closing door	yes
Door hinges	right, reversible
Insulation	70 mm
Gross/net weight	88/83 kg
EAN-Nr.	9005382237199

¹ Indicates the ambient temperature at which the device can be operated safely.



Subject to modification. Timestamp: 09.04.2021



Distributed by:
AMA Medical Products
amamedicalproducts.com.au
sales@amamedicalproducts.com.au
1800 626 292

LIEBHERR