

# LGPv 6520

## MediLine



Gross/net capacity	602/430 litre
Ext. dimensions in mm (w/d/h)	700/830/2160
Int. dimensions in mm (w/d/h)	533/650/1550
Energy consumption in 365 days	1367 kWh
Ambient temperature range	+10°C to +40°C <sup>1</sup>
Refrigerant	R 290
Connection Rating	4.0 A
Voltage	220-240 V~
Cooling system	forced-air
Defrosting	automatic
Temperature range	-9°C to -35°C
Housing material/colour	steel/white
Door material	Steel
Liner Material	304-grade stainless steel
Type of Control	electronic controller
Temperature display	external digital
Fault: warning signal	visual and audible
Interface	RS 485
volt-free alarm contact	yes
Adjustable Shelving	4
Shelves, freezer compartment	4
Castors	castors with brake at the front, fixed castors at the rear
Lock	fitted
Self-closing door	yes
Door hinges	right, reversible
Insulation	83-60 mm
Defrost water drain	yes
Gross/net weight	158/133 kg
EAN-Nr.	9005382197219

<sup>1</sup> Indicates the ambient temperature at which the device can be operated safely.

