

LGv 5010

MediLine



Gross/net capacity	478/337 litre
Ext. dimensions in mm (w/d/h)	751/750/1844
Int. dimensions in mm (w/d/h)	602/520/1224
Energy consumption in 365 days	1245 kWh
Ambient temperature range	+16°C to +35°C ¹
Refrigerant	R 290
Connection Rating	3.0 A
Voltage	220-240 V~
Cooling system	forced-air
Defrosting	automatic
Temperature range	-9°C to -35°C
Housing material/colour	steel/white
Door material	Steel
Liner Material	moulded polystyrol white
Type of Control	electronic controller
Temperature display	external digital
Fault: warning signal	visual and audible
Interface	RS 485
volt-free alarm contact	yes
Drawers	8
Number of baskets	2
Shelves, freezer compartment	5
Castors	castors with brake at the front, fixed castors at the rear
Lock	fitted
Self-closing door	yes
Door hinges	right, reversible
Insulation	70 mm
Defrost water drain	yes
Gross/net weight	124/114 kg
EAN-Nr.	9005382237397

¹ Indicates the ambient temperature at which the device can be operated safely.



Subject to modification. Timestamp: 09.04.2021



Distributed by:
AMA Medical Products
amamedicalproducts.com.au
sales@amamedicalproducts.com.au
1800 626 292

LIEBHERR