

# LKexv 3600

MediLine



Gross/net capacity	333/307 litre
Ext. dimensions in mm (w/d/h)	600/610/1641
Int. dimensions in mm (w/d/h)	470/440/1452
Energy consumption in 365 days	346 kWh
Ambient temperature range	+10°C to +43°C <sup>1</sup>
Refrigerant	R 600a
Connection Rating	1.5 A
Voltage	220-240 V~
Cooling system	forced-air
Defrosting	automatic
Temperature range	+1°C to +15°C
Housing material/colour	steel/white
Door material	Steel
Liner Material	moulded polystyrol white
Type of Control	analogue controls
Temperature display	external digital
Adjustable Shelving	5
Lock	fitted
Door hinges	right, reversible
Insulation	32-52 mm
Gross/net weight	64/59 kg
EAN-Nr.	9005382237496

<sup>1</sup> Indicates the ambient temperature at which the device can be operated safely.



Subject to modification. Timestamp: 09.04.2021



Distributed by:  
AMA Medical Products  
amamedicalproducts.com.au  
sales@amamedicalproducts.com.au  
1800 626 292

# LIEBHERR