

LKexv 5400

MediLine



Gross/net capacity	554/520 litre
Ext. dimensions in mm (w/d/h)	750/729/1640
Int. dimensions in mm (w/d/h)	600/560/1452
Energy consumption in 365 days	359 kWh
Ambient temperature range	+10°C to +40°C ¹
Refrigerant	R 600a
Connection Rating	1.5 A
Voltage	220-240 V~
Cooling system	forced-air
Defrosting	automatic
Temperature range	+1°C to +15°C
Housing material/colour	steel/white
Door material	Steel
Liner Material	moulded polystyrol white
Type of Control	analogue controls
Temperature display	external digital
Adjustable Shelving	5
Lock	fitted
Door hinges	right, reversible
Insulation	50 mm
Gross/net weight	84/77 kg
EAN-Nr.	9005382173992

¹ Indicates the ambient temperature at which the device can be operated safely.

