

# LKv 5710

## MediLine



Gross/net capacity	583/540 litre
Ext. dimensions in mm (w/d/h)	751/750/1844
Int. dimensions in mm (w/d/h)	634/538/1500
Energy consumption in 365 days	438 kWh
Ambient temperature range	+10°C to +35°C <sup>1</sup>
Refrigerant	R 600a
Connection Rating	1.5 A
Voltage	220-240 V~
Cooling system	forced-air
Defrosting	automatic
Temperature range	+3°C to +16°C
Housing material/colour	steel/white
Door material	Steel
Liner Material	moulded polystyrol white
Type of Control	electronic controller
Temperature display	external digital
Fault: warning signal	visual and audible
Interface	RS 485
volt-free alarm contact	yes
Adjustable Shelving	5
Castors	castors with brake at the front, fixed castors at the rear
Lock	fitted
Self-closing door	yes
Door hinges	right, reversible
Insulation	40 mm
Gross/net weight	100/91 kg
EAN-Nr.	9005382238035

<sup>1</sup> Indicates the ambient temperature at which the device can be operated safely.

