

SUFsg 5001

MediLine



Gross/net capacity	491/477 litre
Ext. dimensions in mm (w/d/h)	920/1082/1966
Int. dimensions in mm (w/d/h)	606/605/1300
Energy consumption in 365 days	2884 kWh
Ambient temperature range	+16°C to +32°C ¹
Refrigerant	R 290 / R 290 / R 170
Connection Rating	7.0 A
Voltage	230 V~
Cooling system	static
Defrosting	manual
Temperature range	-40°C to -86°C
Housing material/colour	steel/white
Door material	Steel
Liner Material	stainless steel
Type of Control	electronic controller
Temperature display	external digital
Fault: warning signal	visual and audible
Interface	LAN (Ethernet)
volt-free alarm contact	yes
Adjustable Shelving	3
Shelves, freezer compartment	4
Shelf loading	50 kg
Lock	fitted
Door hinges	right
Insulation	103-180 mm
Defrost water drain	yes
Gross/net weight	308/254 kg
EAN-Nr.	9005382236550

Distributed by:
AMA Medical Products
amamedicalproducts.com.au
sales@amamedicalproducts.com.au
1800 626 292

¹ Indicates the ambient temperature at which the device can be operated safely.

