

# SUFsg 7001

MediLine



|                                      |                             |
|--------------------------------------|-----------------------------|
| <b>Gross/net capacity</b>            | 728/700 litre               |
| <b>Ext. dimensions in mm (w/d/h)</b> | 1204/1082/1966              |
| <b>Int. dimensions in mm (w/d/h)</b> | 890/605/1300                |
| Energy consumption in 365 days       | 2957 kWh                    |
| Ambient temperature range            | +16°C to +32°C <sup>1</sup> |
| Refrigerant                          | R 290 / R 290 / R 170       |
| Connection Rating                    | 7.0 A                       |
| Voltage                              | 230 V~                      |
| Cooling system                       | static                      |
| Defrosting                           | manual                      |
| Temperature range                    | -40°C to -86°C              |
| Housing material/colour              | steel/white                 |
| Door material                        | Steel                       |
| Liner Material                       | stainless steel             |
| Type of Control                      | electronic controller       |
| Temperature display                  | external digital            |
| Fault: warning signal                | visual and audible          |
| Interface                            | LAN (Ethernet)              |
| volt-free alarm contact              | yes                         |
| Adjustable Shelving                  | 3                           |
| Shelves, freezer compartment         | 4                           |
| Shelf loading                        | 50 kg                       |
| Lock                                 | fitted                      |
| Door hinges                          | right                       |
| Insulation                           | 103-180 mm                  |
| Defrost water drain                  | yes                         |
| Gross/net weight                     | 361/313 kg                  |
| EAN-Nr.                              | 9005382236772               |

Distributed by:  
AMA Medical Products  
amamedicalproducts.com.au  
sales@amamedicalproducts.com.au  
1800 626 292

<sup>1</sup> Indicates the ambient temperature at which the device can be operated safely.

