

Chair scales from seca, competence that sits well.

How can you weigh immobile patients or those who are unsteady on their feet as comfortably as possible? It is best to use seca chair scales because they can be brought directly to the patient when needed, sparing unnecessary effort. The intelligent construction with fold-away footrests enables the patient to sit down comfortably and stand up safely. The precise weighing technology provides accurate results.



The practical positioning of the display allows easy operation.

seca 954 Electronic chair scales

Consideration and comfort – these are particularly prominent aspects of this extremely stable model with its circular steel tube frame. A comfortable seat provides a safe position, the swivelling armrests make it easy to move the patient from the bed onto the chair scale and the fold-away footrests increase ease of operation. The locking rear wheels ensure that it is safe to sit down and stand up at all times. The positioning of the large LCD unit behind the chair makes it easy to read and the many functions (including BMI) leave nothing to be desired. This mobile chair scale, with AC-adaptor or rechargeable battery power supply, is an invaluable aid everywhere: in the rehabilitation center, orthopedic department, dialysis center, at nursing homes and in pediatrics and adolescent medicine. For example, small children who are not yet able to stand can be weighed comfortably while sitting. The practical pouch seca 471 attached to the scale stores the mains adapter.

- Capacity: 200 kg
- Graduation: 100 g
- Functions: Automatic switch-off, TARE, HOLD, BMI, kg/lbs/sts switch-over
- Optional: Available with integrated PC or printer interface



To make sitting down and standing up easier, the arm and footrests can be folded away.

