

# seca supra 719

A sculpted form that endures the test of time.



- Display on a high column.  
For easy reading.
- Step-off function. For weighing without switching on scales.
- High stability.  
For demanding users.
- More than 170 years of experience. For maximum precision.
- From the global market leader for medical weighing and measuring.

## seca supra 719:

A sculpted form that endures the test of time.

**T**ake a look at the four load cells. Made from finely brushed stainless steel, they ensure a stable connection between weighing platform and base.



### A sculpted form that radiates precision.

Slim, sleek, simple. Who would believe that this scale can do more than just look good? High-quality ABS plastic (acrylonitrile-butadiene-styrene) withstands impacts, shocks and all commercially available cleaning or disinfecting products. High-grade 8-millimetre-thick safety glass ensures a good feeling. Made for people who work on their bodies or who have to pay attention to their health and therefore rely on precise measurements. It's all possible at home with these scales, the most precise scales for private use.

### A sculpted form that radiates beauty.


The six-screw connection between base and column ensures maximum solidity. After all, it looks good in the sales room and after many years of use at home.

### A sculpted form that radiates simplicity.

What is simpler than the Step-off function? Simply stand on the scale without switching it on, wait for "Step-off" on the display, then Step off at leisure and read off the measurements. This is ideal particularly for pregnant women and people of larger girth. You see, a scale from seca is on the ball. It's just one of the reasons why seca has been the global market leader in the field of medical weighing and measuring for many years.



*The display on a high column is easy to read.*



*Four load cells made from brushed stainless steel.*

## seca supra 719

### Technical data

- Capacity: 180 kg / 400 lbs
- Division: 100 g / 0.2 lbs
- Dimensions (WxHxD):  
350 x 955 x 435 mm / 13.8 x 37.6 x 17.1"
- Weight: 6.1 kg / 13.4 lbs
- Functions: Step-off, auto-HOLD, automatic switch-off, tip-on function, kg/lbs switch-over
- Parts to be delivered: 4 AAA batteries, 4 extension feet for carpets

