

seca integration solutions
A clear improvement
to your workflow

The seca 787 and seca 797 can be integrated into your information system via the existing network. At the push of a button, measured values can be sent directly to the electronic medical record (EMR), where they are immediately available to all users.

Connected with an optional USB barcode scanner, both patient and user are identified and the measurement results are always correctly assigned. This saves time, increases process reliability and eliminates transcription errors.

Get further information on seca integration solutions at:
integration.seca.com

The seca integration workflow is supported by several system providers:



seca 777 series

seca

Precision for health

Time to say good buy!



Technical data

	seca 777	seca 787	seca 797
Functions	TARE, HOLD, lbs/kg switch-over lock, BMI	TARE, HOLD, Auto-HOLD, Damping (motion filter), lbs/kg switch-over lock, automatic BMI, auto-off	TARE, HOLD, Auto-HOLD, Damping (motion filter), kg/lbs switch-over lock, automatic BMI, USB-port and holder for scanner
Data transmission	USB interface for customer specific solution	LAN or Wi-Fi via interface module seca 452 (optional)	Integrated Wi-Fi module
Power supply	Batteries, power adapter (optional)	Batteries, power adapter (optional)	Power adapter, battery pack (optional)
Weight	12.7 kg / 28.0 lbs	12.9 kg / 28.4 lbs	13.1 kg / 28.9 lbs
Measuring range (total)	10 – 230 cm / 4 – 90 inch	10 – 230 cm / 4 – 90 inch	10 – 230 cm / 4 – 90 inch
Measuring range (digital)		122 – 230 cm / 48 – 90 inch	122 – 230 cm / 48 – 90 inch
Capacity		250 kg / 550 lbs	
Graduation		100 g / 0.2 lbs	
Dimensions (W x H x D)		41.0 x 135.6 x 57.6 cm / 16.1 x 53.4 x 22.7 inch	
Platform (W x H x D)		35.0 x 5.9 x 41.8 cm / 13.8 x 2.3 x 16.5 inch	
Electromagnetic radiation resistance		60601-1-2 compliant	

seca headquarters
+49 40 20 00 00 0
info@seca.com
Hamburg, Germany

seca middle east
+971 564 410 292
info.ae@seca.com
Dubai, United Arab Emirates

seca asia pacific
+60 173 711 973
info.aspac@seca.com
Kuala Lumpur, Malaysia

Precision for health
seca.com

