

MEDICAL PRODUCTS  
**WALMAY**

# CRYOSPRAY

MODEL - WM161

The Walmay Cryospray has given over two decades of well proven reliable service to Medical Professionals.

It is the instrument of choice for all Dermatologists and Practitioners.

WALMAY CRYOSPRAY



NEW & IMPROVED MODEL

FIVE YEAR WARRANTY

Proven features that makes the Walmay Cryospray unit WM161 the preferred choice for Medical Professionals around the world.

**5 YEAR  
WARRANTY**

Fixed Stainless Steel Probe with greater tensile strength.

Increased travel in the Regulator Handle allowing better control with the flow of Liquid Nitrogen

The centre of gravity of the Walmay Cryospray WM161 is towards the top of the flask, making it easy for the doctor to balance and steady the instrument in operation.

Includes sturdy protective Carry Bag containing five brass spray tips, user and maintenance instructions, and FIVE YEAR Warranty card.

Improved metallurgical principles in construction of the flask collar.

Fine gauze Strainer on the end of the Siphon Tube to reduce contamination accessing the Regulator mechanism.

The base stand provides a greater surface area for stability of the instrument.

Includes USB containing a video for setting up and using the Walmay Cryospray WM161.



AUSTRALIAN  
DESIGN &  
MANUFACTURE

ARTG  
NUMBER  
170073

EUROPEAN  
APPROVAL  
CE 0805