

# seca 664

## Digital wheelchair scale with wireless transmission

**EMR integrated**



- Network-capable with seca 360° wireless technology.
- High capacity, up to 800 lbs.
- Flat platform with self-lowering ramp.
- Two fold-away ramps included.
- Easy-to-fold and transport.



Distributed by:  
AMA Medical Products  
amamedicalproducts.com.au  
sales@amamedicalproducts.com.au  
1800 626 292

**seca**<sup>®</sup>  
Precision for health

## seca 664:

Real comfort with maximum safety.

**S**afety, functionality and aesthetic design – the seca 664 sets the standard in medical weighing. Specially designed for weighing wheelchair patients in rehabilitation centers, nursing homes, hospitals and in dialysis centers – this scale provides both its users and patients maximum comfort and safety.

### High capacity and easy to mount ramp.

With the scale's load-bearing capacity of 800 lbs, weighing even heavy wheelchair patients is effortless with the seca 664. An extremely flat platform with an automatically lowered ramp ensure that the scale is easy to mount. The firmly installed side rails and a black ribbed rubber surface are additional safety features.

### Integrated pre-TARE function.

With the pre-TARE function, the weight of the wheelchair can first be stored. Then, at the touch of a button, the weight of the wheelchair can be deducted from the total weight of the patient and wheelchair, thus automatically calculating the weight of the patient.

### Paperless and error-free wireless transmission of measurements.

With the seca 360° wireless technology, the seca 664 can transmit measurements to the optional seca 360° wireless digital printer. With the network-capable seca analytics software and seca 456 USB adapter, your PC can receive and analyze measurements. The seca 664 is EMR-integrated and ready to handle all the demands the future will bring. For more information on EMR integration, visit: [www.seca.com](http://www.seca.com)



## seca 664

### Technical Data

- Capacity: 800 lbs / 360 kg
- Graduation: 0.1 lbs / 50 g
- Dimensions (WxHxD): 43.4 x 35.9 x 45.3" / 1,102 x 912 x 1,150 mm
- Dimensions/platform (WxHxD): 31.5 x 2.2 x 38" / 800 x 55 x 965 mm
- Weight: 73.6 lbs / 33.4 kg
- Power supply: power adapter, rechargeable batteries
- Functions: TARE, pre-TARE, HOLD, auto-HOLD, BMI, auto-CLEAR, lbs/kg switch-over, damping, automatic switch-off, SEND/PRINT
- Included: Two fold-away ramps
- RS232 compatible
- System compatible with: 360° wireless digital printer advanced seca 466, 360° wireless digital printer seca 465, 360° wireless USB adapter seca 456, software seca analytics 115 and seca emr flash 101, stadiometer seca 274 or seca 264



Distributed by:  
AMA Medical Products  
[amamedicalproducts.com.au](http://amamedicalproducts.com.au)  
[sales@amamedicalproducts.com.au](mailto:sales@amamedicalproducts.com.au)  
1800 626 292